Spesera

Health Products, Imp. & Exp. Trade Ltd.

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Ref. No.: MU-EAR-22101

Date: 09/03/2022

This is to certify that, according to the Council Directive 98/79/EC, MedUnion S.L performed all notification duties and responsibilities as the European Authorized Representatives (EC REP) of:

Manufacturer: Spesera Health Products. Imp & Exp. Trade, Ltd
Address: Goztepe mah. Batisehir cad. No:2/3 Bagcilar/ISTANBUL

The Manufacturer has provided MedUnion S.L. with all the appropriate declarations according to the 98/79/EC Directive-article 10 requirements including the EC Declaration of Conformity confirming that his In-Vitro Diagnostics medical devices, as stipulated in the Annex I, are fulfilling the applicable requirements of the European Council Directive 98/79/EC

The notification of the In-Vitro Diagnostic medical devices in the Annex I has been completed by MedUnion S.L on the **09/03/2022** with Registratation number **RPS/602/2022** in Spanish Medicines and Medical Products Agency.

Conformity Assessment Route: Medical device directive: Council Directive 98/79/EC concerning in vitro diagnostic medical devices (IVDD 98/79/EC).

In-Vitro Diagnostic medical devices: Please see ANNEX I-List of devices (1 page, 22 devices)

Ms Shu Chuan Kuo, Managing Director
MedUnion S.L.

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
BELLA Certification certify that the Quality Management Systems of the above organization has been audited and found to be in accordance with the requirements of the management system standards detailed below Standard

ISO 9001:2015

Scope of Certification

Packaging of Mikrobiological Products, Preparation of Medical Laboratory Materials,
Reception, Sale, Import and Export

EA CODE: 29

Certificate No : 19-217-1115
Certification Cycle Start Date : 26.03.2019
Reissue Date : 21.03.2022/03
Certificate Expire Date : 25.03.2025

Validity Date of Current Certification Period: 21.03.2022 – 25.03.2023

TÜRKAK BDS NO
YS-A685-5D57

This certificate is valid for above address(es) and scope upon compliance with Bella procedure(s), instruction(s) and the above mentioned standard and successful completion of regular surveillance audits. The validity of this certificate can be checked through on www.bqcert.co.uk webpage.
Bella Kalite Belgelendirme Hiz. Ltd. Şti. Cevizli Mah. Zuhul Cad. No:44 B Blok D:56
Ritim İstanbul Konutları Maltepe/İSTANBUL - Tel: 0 216 518 33 08 Fax: 0 216 518 40 74 www.bqcert.co.uk - info@bqcert.co.uk

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BELLA has evaluated and approved the suitability of the Microbiology of Food, Animal Feed and Water Preparation, Production, Storage and Performance Testing of Culter Media Management Systems within the scope of the following.

ISO 11133:2014

Production, Import, Marketing and Export of Devices and Chemicals
Equipments That Used in the Fields of Molecular Biology,
Microbiology and Biotechnology

Certificate No : 19-19005
Certification Cycle Start Date : 18.03.2019
Reissue Date : 18.03.2022/03
Certificate Expire Date : 18.03.2023
Validity Date : 1 Year

This certificate is valid for above address(es) and scope upon compliance with Bella procedure(s), instruction(s) and the above mentioned standard and successful completion of regular surveillance audits. The validity of this certificate can be checked through on www.bqcert.co.uk webpage.
Bella Kalite Belgelendirme Hiz. Ltd. Şti. Cevizli Mah. Zuhai Cad. No:44 B Blok D:56
Ritim Istanbul Konutları Maltepe/İSTANBUL - Tel: 0 216 518 33 08 Fax:0 216 518 40 74 www.bqcert.co.uk - info@bqcert.co.uk

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Göztepe Mah. Batişehir Cad. Numara:2/3 No:1B 12 Bağcılar /İstanbul

BELLA has evaluated and approved the compliance of the Medical Devices
Management Systems within the scope of the following

ISO 13485:2016

Production, Import, Marketing and Export of Devices and Chemicals
Equipments That Used in the Fields of Molecular Biology,
Microbiology and Biotechnology

Certificate No : 19-19006
Certification Cycle Start Date : 18.03.2019
Reissue Date : 18.03.2022/03
Certificate Expire Date : 18.03.2023
Validity Date : 1 Year

This certificate is valid for above address(es) and scope upon compliance with Bella procedure(s), instruction(s)
and the above mentioned standard and successful completion of regular surveillance audits. The validity of this certificate
can be checked through on www.bqcert.co.uk webpage.
Bella Kalite Belgelendirme Hiz. Ltd. Şti. Cevizli Mah. Zuhul Cad. No:44 B Blok D:56
Ritim İstanbul Konutları Maltepe/İSTANBUL - Tel: 0 216 518 33 08 Fax:0 216 518 40 74 www.bqcert.co.uk - info@bqcert.co.uk

FR.12/03.02.2021/04

**THE REPUBLIC OF TURKEY
MINISTRY OF HEALTH
TURKISH MEDICINES AND MEDICAL DEVICES AGENCY**

Certificate No: 163897

Date of Issue : 17 September 2021

CERTIFICATE OF FREE SALE

To whom it may concern,

It is hereby certified that the products detailed in the attached schedule, which are manufactured by "SPESERA SAĞLIK ÜRÜNLERİ İÇ VE DIŞ TİC. LTD. ŞTİ." (Göztepe Mh. Batışehir Cd. No:2/3/1b12 BAĞCILAR İSTANBUL), have been affixed with the CE mark in accordance with Medical Device Directives of the European Union (EU) and are freely sold in Turkey.

This certificate is issued to be given to the relevant competent authorities of "Bulgaria" and is valid for 36 months from the date of issue.

Yours sincerely,

Ömer Faruk KURU
Head of Medical Devices
Registration and Coordination Department

This certificate consists of 4 page/s and 22 products. The products listed in the attached schedule are registered from the date of issuance of this certificate and information about the current status of these products is accessible through

<https://utsuygulama.saglik.gov.tr/UTS/vatandas#/vatTibbiCihazListele>.

PRODUCT SCHEDULE

#	Barkod	Brand	Label Name	Reference No / Version / Model	GMDN Code
1	8681176311002	SPESERA	COVID-19 IgG/IgM Rapid Test Cassette	SP-225	64756
2	8681176311019	SPESERA	sNAT Transfer Tube	SP-227	58697
3	8681176311026	SPESERA	SPESERA COVID - 19 Antigen Rapid Test Cassette	SP-226	64787
4	8681176311125	SPESERA	Influenza A+B Rapid Test Cassette	SP-222	49119
5	8681176311156	SPESERA	H.PYLORI ANTIGEN Rapid Test Cassette	SP-207	30825
6	8681176311354	SPESERA	FOB Rapid Test Cassette	SP-304	54532
7	8681176311385	SPESERA	HAV IgG/IgM Rapid Test Cassette	SP-215	48270
8	8681176311064	SPESERA	TROPONIN I Rapid Test Cassette	SP-306	46989
9	8681176311095	SPESERA	CHLAMYDIA Rapid Test Cassette	SP-109	50759
10	8681176311101	SPESERA	CK-MB Rapid Test Cassette	SP-307	52995

11	8681176311088	SPESERA	BETA HCG (PREGNANCY) Rapid Test Cassette	SP-103	33819
12	8681176311132	SPESERA	Beta HCG DIPSTICK	SP-101	33819
13	8681176311149	SPESERA	Urine Strip 10T	SP-309	54514
14	8681176311187	SPESERA	SYPHILIS Rapid Test Cassette	SP-219	51788
15	8681176311194	SPESERA	GIARDIA LAMBLIA Rapid Test Cassette	SP-205	52249
16	8681176311248	SPESERA	H.PYLORI ANTIBODY Rapid Test Cassette	SP-206	62029
17	8681176311231	SPESERA	STREP A HIZLI TEST	SP-220	51707
18	8681176311224	SPESERA	DENGUE IgG/IgM Rapid Test Cassette	SP-208	63238
19	8681176311200	SPESERA	TYPHOID Rapid Test Cassette	SP-211	51560
20	8681176311347	SPESERA	TUBERCULOSIS (TB) Rapid Test Cassette	SP-223	51172
21	8681176311279	SPESERA	ADENOVİRÜS Rapid Test Cassette	SP-201	49856

Date of Issue : 17 September 2021

22	8681176311170	SPESERA	ROTAVİRÜS- ADENOVİRÜS Rapid Test Cassette	SP-204	48235
----	---------------	---------	--	--------	-------

End of product schedule.

Live, work, travel in the EU

COVID-19 In Vitro Diagnostic Devices and Test Methods Database

[Home](#) > [COVID-19 In Vitro Diagnostic Medical Devices](#) >

COVID-19 In Vitro Diagnostic Medical Device - detail

COVID-19 In Vitro Diagnostic Medical Device - detail

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Device identification number	2246
CE Marking	Yes
HSC common list	No
Format	Near POC / POC
Physical Support	Cassette
Target type	Antigen
Specimen	Nasal swab, Nasopharyngeal swab, Throat swab

Cross-reactivity (pathogens tested)	Alpha Coronavirus 229E (HCoV-229E), Alpha Coronavirus NL63 (HCoV-NL63), Beta Coronavirus HKU1 (HCoV-HKU1), Beta Coronavirus OC43 (HCoV-OC43)
Commercial Status	Commercialised
Last Update	2022-03-14 09:42:40 CET
Comments	Covid-19 Antigen Rapid Test Cassette

Assay Type	Immuno-Antigen
Reader Required	No
Subcategory	Extraction kit
Method	Immunochromatography
Measurement	Qualitative
Time	15 minutes
Subclass	Multiple
LOD	580 TCID50/ml
Positive control	$100 \times TP / (TP + FP) = 100 \times 238 / 242 = 98.34\%$
Negative control	$100 \times TN / (FN + TN) = 100 \times 530 / 538 = 98.51\%$
Analytical Sensitivity	93.76 %
Analitical Specificity	97.5 %
Clinical Sensitivity	96.89 %
Clinical Specificity	99.1 %
Notes	Ares(2021)3875860

The database contains publicly available In Vitro Diagnostic Medical Devices for COVID-19 and it is being updated periodically. Please note that additional performance

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URINALYSIS REAGENT STRIPS


Health Products, Imp. & Exp. Trade Ltd.

Drugs of Abuse/Tumor & Cardiac Markers /Other PRODUCT LIST

Cat. No.	Product Description	Specimen	Format	Kit Size	Status
SP-301	Multi-drug 6drugs Test Cassette	Urine	Cassette	25 T	CE
SP-302	Multi-drug 10drugs Test Cassette	Urine	Cassette	25T	CE
SP-303	Multi-drug 12 drugs Test Cassette	Urine	Cassette	25 T	CE
SP-304	FOB Rapid Test Cassette	Feces	Cassette	25 T	CE
SP-305	PSA Rapid Test Cassette	WB/S/P	Cassette	25 T	ISO
SP-306	cTnI Rapid Test Cassette	WB/S/P	Cassette	25 T	CE
SP-307	CK-MB Rapid Test Cassette	WB/S/P	Cassette	10 T/25 T	CE
SP-308	Vitamin D Rapid Test Cassette	WB	Cassette	10 T	CE
SP-309	Urinalysis Reagent Dipsticks - 10 parameters	Urine	Dipstick	100 T	CE


sNAT TRANSFER TUBE FOR COVID-19

TRANSFER TUBE FOR COVID-19

Nasopharyngeal Swab + transfer Tube

NUCLISWAB* is a sampling and transport system used for collection and transport of viruses and cells from clinical samples or from environment, for isolation of nucleic acids that will be used in nucleic acid amplification tests.

Spesera
Health Products, Imp. & Exp. Trade Ltd.

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
TEST PROCEDURE

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COVID-19 IgG/IgM Rapid Test

COVID-19 IgG/IgM Rapid Test Cassette (Whole Blood/Serum/Plasma) is a solid phase immunochromatographic assay for the rapid, qualitative and differential detection of IgG and IgM antibodies to 2019 Novel Coronavirus in human whole blood, serum or plasma.

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Tel: 00905350499266
www.spesera.com.tr


Key Facts of COVID-19

- The 2019 new coronavirus pneumonia (NCP), or “COVID-19”, was discovered for the 2019 Wuhan Viral Pneumonia case in China and was named by World Health Organization on January 12, 2020.
- For confirmed coronavirus disease 2019 (COVID-19) cases, reported illnesses have ranged from mild symptoms to severe illness and death. Symptoms can include:
 - Fever
 - Cough
 - Shortness of breath
- The virus is thought to spread mainly from person-to-person.
 - Between people who are in close contact with one another (within about 6 feet).
 - Via respiratory droplets produced when an infected person coughs or sneezes.
 - These droplets can land in the mouths or noses of people who are nearby or possibly be inhaled into the lungs.

Features & Benefits

- Fast results as soon as 2-10 minutes.
- Facilitates patient treatment decisions quickly.
- Simple, time-saving procedure.
- Little specimens, only 5 µL of serum/plasma or 10 µL of whole blood specimens.
- All necessary reagents provided & no equipment needed.
- High sensitivity and specificity.

Test Procedure & Interpretation

Ordering Information

Product Description	Specimen	Catalog No.	Format	Kit Size
COVID-19 IgG/IgM Rapid Test	Whole Blood/Serum/Plasma	GCCOV-402a	Cassette	25 Tests/Kit

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Health Products, Imp. & Exp. Trade Ltd.

COVID-19 ANTIGEN RAPID TEST KIT (ORAL FLUID)

Detection kit for SARS-CoV-2 antigen in sputum or Posterior oropharyngeal saliva from human


INTENDED USE

COVID-19 Antigen Rapid Test Kit is an immunochromatographic assay kit for the qualitative detection of SARS-CoV -2 antigen in sputum and Posterior oropharyngeal saliva from human.

TEST PROCEDURE

Please read the instructions carefully and allow the test device and specimens to equilibrate to temperature (15°C-30°C) prior to testing. Before performing the test, please read the instructions carefully and allow the test device and sample to reach the indoor temperature (15°C-30°C)

1. Saliva and sputum specimens from the back of the oropharynx: Use soap and water mixed with alcohol to clean the hands. The best time for specimen collection: after getting up, before brushing your teeth, eating or drinking.
2. Take out the saliva sampling bag, make a "Kua Kua" sound from the deep throat, clear the saliva from the deep throat, and then spit the saliva (about 2ml) into the bag. Do not have any saliva contamination on the surface of the sampling bag.
3. Use the dropper provided to stir the sputum sample 8-10 times, then aspirate the sputum sample and add it to the extraction tube.
4. Take out the extraction tube containing the extraction buffer then mix the extraction buffer and specimen well.
5. Take out the test device from the packaging box and place it on a flat table.
6. Add 4 drops of the mixed solution to the sample hole of the test card.
7. Read the result at 15 minutes. If you haven't read it for 20 minutes or more, the result is invalid and it is recommended to repeat the test..

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

Spesera

Health prouducts, Imp. & Exp. Trade Ltd.

Influenza & COVID-19 Ag Combo Rapid Test Cassette (Swab)

INTENDED USE

The Influenza & COVID-19 Ag Combo Rapid Test Cassette (Swab) is an in vitro immunochromatographic assay for the qualitative and differential detection of nucleocapsid protein antigen from influenza A (including the subtype H1N1), influenza B and/or SARS-CoV-2 in nasopharyngeal (NP) swab specimens. It is intended to aid in the rapid diagnosis of influenza A, influenza B and/or SARS-CoV-2 infections. This test provides only a preliminary test result. Therefore, any reactive specimen with the Influenza & COVID-19 Ag Combo Rapid Test Cassette (Swab) must be confirmed with alternative testing method(s) and clinical findings.

SUMMARY AND EXPLANATION

Influenza is an acute and highly contagious viral infection of the respiratory tract. The causative agents of the disease are immunologically diverse, single-strand RNA virus known as influenza viruses. There are three types of influenza viruses: A, B and C. Type A viruses are the most prevalent and are associated with most serious epidemics, while Type B infection is generally milder. Type C virus have never been associated with a large epidemic of human disease. Both type A and B viruses can circulate simultaneously, but usually one type is dominant during a given season and particular epidemic area. The disease is easily transmitted through coughing and sneezing of aerosolized droplets containing live virus. Influenza outbreaks normally occur each year during fall and winter seasons.

The novel coronaviruses belong to the β genus. COVID-19 is an acute respiratory infectious disease. Entire human population is susceptible. Currently, the patients infected by the novel coronavirus are the main source of infection; asymptomatic infected people can also be an infectious source. Based on the current epidemiological investigation, the incubation period is 1 to 14 days, mostly 3 to 7 days. The main manifestations include fever, fatigue, dry cough, and loss of taste and smell. Nasal congestion, runny nose, sore throat, myalgia and diarrhea are found in a few cases. This test is for detection of SARS-CoV-2 nucleocapsid protein antigen. Antigen is generally detectable in upper respiratory specimens during the acute phase of infection. Rapid diagnosis of SARS-CoV-2 infection will help healthcare professionals to treat patients and control the disease more efficiently and effectively.

PRINCIPLE OF THE TEST

The Flu A&B Antigen Rapid Test is an immunochromatographic membrane assay that uses highly sensitive monoclonal antibodies to detect influenza type A and B nucleoprotein antigens in nasopharyngeal swab samples. The test strip is composed of the following parts: namely sample pad, reagent pad, reaction membrane, and absorbing pad. The reagent pad contains the colloidal-gold conjugated with the monoclonal antibodies against Influenza virus A and B; the reaction membrane contains the secondary antibodies either for virus A or for B. The whole strip is fixed inside a plastic device. When the sample is added into the sample well, conjugates dried in the reagent pad are dissolved and migrate along with the sample. If influenza A presents in the sample, a complex formed between the anti-influenza A conjugate and the virus will be captured by the specific anti-influenza A monoclonal antibodies coated on the A region (A). If the sample contains influenza B, a complex formed between the anti-influenza B conjugate and the virus will be captured by the specific anti-influenza B monoclonal antibodies coated on the B region (B). Results appear at 10 minutes in the form of a red line that develops on the membrane. To serve as a procedural control, a red line will always appear in the control region (C) indicating that proper volume of sample has been added and membrane wicking has occurred.

The COVID-19 Antigen Rapid Test Cassette is an immunochromatographic membrane assay that uses highly sensitive monoclonal antibodies to detect nucleocapsid protein from SARS-CoV-2 in nasopharyngeal (NP) swab. The test strip is composed of the following parts: namely sample pad, reagent pad, reaction membrane, and absorbing pad. The reagent pad contains the colloidal-gold conjugated with the monoclonal antibodies against the nucleocapsid protein of SARS-CoV-2; the reaction membrane contains the secondary antibodies for nucleocapsid protein of SARS-CoV-2. The whole strip is fixed inside a plastic device. When the sample is added into the sample well, conjugates dried in the reagent pad are dissolved and migrate along with the sample. If SARS-CoV-2 antigen presents in the sample, a complex formed between the anti-SARS-2 conjugate and the virus will be captured by the specific anti-SARS-2 monoclonal antibodies coated on the test line region (T). Absence of the T line suggests a negative result. To serve as a procedural control, a red line will always appear in the control line region (C) indicating that proper volume of sample has been added

and membrane wicking has occurred.

MATERIALS PROVIDED

1. 20 Test cassettes
2. 20 Sterile swabs
3. 20 Extraction tubes and 20 dropper tips
4. 1 Workstation
5. 2 Buffers
6. 1 Package insert

MATERIALS REQUIRED BUT NOT PROVIDED

- 1.Clock,timer or stopwatch

WARNINGS AND PRECAUTIONS

- 1.For in vitro diagnostic use only.
- 2.The test device should remain in the sealed pouch until use.
- 3.Do not use kit past its expiration date.
- 4.Swabs, tubes and test devices are for single use only.
- 5.The extraction buffer contains a solution with a preservative (0.09% sodium azide). If solution comes in contact with the skin or eyes, flush with ample volumes of water.
- 6.Solutions that contain sodium azide may react explosively with lead or copper plumbing. Use large quantities of water to flush discarded solutions down a sink.
- 7.Do not interchange or mix components from different kit lots.
- 8.When collecting a nasopharyngeal swab sample, use the Nasopharyngeal Swab supplied in the kit.
- 9.Wear appropriate personal protection equipment and gloves when running each test and handling patient specimens. Change gloves between handling of specimens suspected of COVID-19.
- 10.Specimens must be processed as indicated in the SPECIMEN COLLECTION and SAMPLE PREPARATION

PROCEDURE sections of this Product Insert. Failure to follow the instructions for

use can result in inaccurate results.

- 11.To obtain accurate results, do not use visually bloody or overly viscous samples.
- 12.Proper laboratory safety techniques should be followed at all times when working with SARS-CoV-2 patient samples. Patient swabs, used Test Strips and used extraction buffer vials may be potentially infectious. Proper handling and disposal methods should be established by the laboratory in accordance with local regulatory requirements.
- 13.Humidity and temperature can adversely affect results.
- 14.Used testing materials should be discarded in accordance with local regulations.

STORAGE AND STABILITY

1. The kit can be stored at room temperature or refrigerated (2-30°C).
2. Do not freeze any of the test kit components.
3. Do not use test device and reagents after expiration date.
4. Test devices that have been outside of the sealed pouch for more than 1 hour should be discarded.

SPECIMEN COLLECTION


Use the nasopharyngeal swab supplied in the kit.

- 1.Carefully insert the swab into the nostril of the patient, reaching the surface of posterior nasopharynx. that presents the most secretion under visual inspection.
- 2.Swab over the surface of the posterior nasopharynx. Rotate the swab several times.
- 3.Withdraw the swab from the nasal cavity.


SAMPLE PREPARATION PROCEDURE

1. Insert the test extraction tube into the paper stand in this product. Make sure that the tube is standing firm and reaches the bottom of the stand.
2. Add the sample extraction buffer to the extraction tube until it reaches the lower mark (about 13-17 drops, 0.5 mL).

3. Insert the swab into the extraction tube which contains 0.5 ml of the extraction buffer.
4. Roll the swab at least 6 times while pressing the head against the bottom and side of the extraction tube.
5. Leave the swab in the extraction tube for 1 minute.
6. Squeeze the tube several times with fingers from outside of the tube to immerse the swab. Remove the swab. The extracted solution will be used as test sample.

SPECIMEN TRANSPORT AND STORAGE

Do not return the nasopharyngeal swab to the original paper packaging.


For best performance, direct nasopharyngeal swabs should be tested as soon as possible after collection. If immediate testing is not possible, and to maintain best performance and avoid possible contamination, it is highly recommended the nasopharyngeal swab is placed in a clean, unused plastic tube labeled with patient information, preserving sample integrity, and capped tightly at room temperature

(15-30°C) for up to 1 hour prior to testing. Ensure the swab fits securely within the tube and the cap is tightly closed. If greater than 1 hour delay occurs, dispose of sample. A new sample must be collected for testing.

TEST PROCEDURE

Allow the test device, test sample and buffer to equilibrate to room temperature (15-30°C) prior to testing.

1. Remove test device from the sealed pouch just prior to the testing and lay flat on work bench.
2. Insert a nozzle with filter into the sample extraction tube tightly.
3. Reverse the sample extraction tube, and add 4 drops (about 100 µl) of test sample by squeezing the extracted solution tube into the each sample well (S).
4. Wait for the colored band(s) to appear. The result should be read in 15 minutes. Do not interpret the result after 20 minutes.

INTERPRETATION OF RESULTS

For Flu A&B Antigen Rapid Test

1. POSITIVE:

1.1 Flu A Positive:

The presence of two lines as control line (C) and A test line within the result window indicates a positive result for Influenza A viral antigen.

1.2 Flu B Positive:

The presence of two lines as control line (C) and B test line within the result window indicates a positive result for Influenza B viral antigen.

1.3 Flu A+B Positive:

The presence of three lines as control line (C), A test line and B test line within the result window indicates a positive result for Influenza A and Influenza B viral antigen.

2. NEGATIVE:

The presence of only control band (C) within the result window indicates a negative result.

3. INVALID:

If the control band (C) is not visible within the result window after performing the test, the result is considered invalid. Some causes of invalid results are because of not following the directions correctly or the test may have deteriorated beyond the expiration date. It is recommended that the specimen be re-tested using a new test.

For COVID-19 Antigen Rapid Test

1. POSITIVE:

The presence of two lines as control line (C) and test line (T) within the result window indicates a positive result.

2. NEGATIVE:

The presence of only control line (C) within the result window indicates a negative result.

3. INVALID:

If the control line (C) is not visible within the result window after performing the test, the result is considered invalid. Some causes of invalid results are because of not following the directions correctly or the test may have deteriorated beyond the expiration date. It is recommended that the specimen be re-tested using a new test.

NOTE:

1.The intensity of color in the test line region (T) may vary depending on the concentration of analyses present in the specimen. Therefore, any shade of color in the test line region (T) should be considered positive. Please note that this is a qualitative test only, and cannot determine the concentration of analytes in the specimen.

2.Insufficient specimen volume, incorrect operating procedure or expired tests are

the most likely reasons for control band failure.

QUALITY CONTROL

A procedural control is included in the test. A red line appearing in the control line region (C) is the internal procedural control. It confirms sufficient specimen volume and correct procedural technique. Control standards are not supplied with this test. However, it is recommended that positive and negative controls are sourced from a local competent authority and tested as a good laboratory practice, to confirm the test procedure and verify the test performance.

LIMITATIONS

1. The Influenza & COVID-19 Ag Combo Rapid Test Cassette (Swab) is for professional in vitro diagnostic use, and should only be used for the qualitative detection of influenza A ,B and/or SARS-CoV-2 in nasopharyngeal (NP) swab specimens.

2. The etiology of respiratory infection caused by microorganisms other than influenza A , B or SARS-CoV-2 will not be established with this test.

3. The Influenza & COVID-19 Ag Combo Rapid Test Cassette (Swab) is capable of detecting both viable and non-viable influenza and SARS-CoV-2 particles. The performance of the Influenza & COVID-19 Ag Combo Rapid Test Cassette (Swab) depends on antigen load and may not correlate with cell culture performed on the same specimen.

4. If the test result is negative and clinical symptoms persist, additional testing using other clinical methods is recommended. A negative result does not at anytime rule out the presence of influenza A A , B or SARS-CoV-2 viral antigens in specimen, as they may be present below the minimum detection level of the test. As with all diagnostic tests, a confirmed diagnosis should only be made by a physician after all clinical and laboratory findings have been evaluated.

5.Inadequate or inappropriate specimen collection, storage, and transport may yield false negative test result.

6. Failure to follow the Test Procedure may adversely affect test performance and/or invalidate the test result.

7. Although this test has been shown to detect cultured avian influenza viruses, including avian influenza A subtype H5N1 virus, the performance characteristics of

this test with specimens from humans infected with H5N1 or other avian influenza viruses are unknown.

8. Performance characteristics for influenza A were established when influenza A/H3 and A/H1 were the predominant influenza A viruses in circulation. When other influenza A viruses are emerging, performance characteristics may vary.

9. Positive and negative predictive values are highly dependent on prevalence. False positive test results are more likely during periods of low influenza activity when prevalence is moderate to low. Positive test results do not rule out co-infections with other pathogens.

10. For COVID-19 Antigen Rapid Test, positive test results do not differentiate between SARS-CoV and SARS-CoV-2. Negative results should be treated as presumptive and confirmed with an FDA authorized molecular assay, if necessary, for clinical management, including infection control.

PERFORMANCE CHARACTERISTICS

For FluA&B Antigen Rapid Test:

1. Analytical Sensitivity

The minimum detection limit is 1.5×10^4 TCID₅₀ /test for the Influenza A virus antigen and is 1.5×10^5 TCID₅₀ /test for the Influenza B virus antigen.

2. Analytical Reactivity

The influenza A strain listed tested positive in the Influenza A & B Ag Rapid Test. Although the specific influenza strains causing infection in human can vary, all contain the conserved nucleoproteins targeted by Influenza A&B Ag Rapid Test.

Strains	Sources	Subtypes	Concentration
Flu A/Hubei/PR8/2001	Human	H1N1	1.8×10^4 TCID ₅₀ /test
Flu A/New Kaledonia/20/99	Human	H1N1	1.8×10^4 TCID ₅₀ /test
Flu A/Yamagata/32/89	Human	H1N1	1.8×10^4 TCID ₅₀ /test
Flu A/Beijing/262/95	Human	H1N1	1.8×10^4 TCID ₅₀ /test
Flu A/Singapore/1/57	Human	H2N2	3.0×10^4 TCID ₅₀ /test
Flu A/Hubei/3/2005	Human	H3N2	3.0×10^4 TCID ₅₀ /test
Flu A/Akita/1/94	Human	H3N2	3.0×10^4 TCID ₅₀ /test
Flu A/Kita Kyusyu/159/93	Human	H3N2	3.0×10^4 TCID ₅₀ /test
Flu A/Iowa/15/30	Swine	H1N1	3.0×10^4 TCID ₅₀ /test
Flu A/Hongkong/168/93	Swine	H1N1	3.0×10^4 TCID ₅₀ /test
Flu A/Anhui/24/2004	Swine	H5N1	6.0×10^4 TCID ₅₀ /test
Flu A/Hubei/134/2000	Swine	H9N2	6.0×10^4 TCID ₅₀ /test
Flu A/Hubei/251/2001	Swine	H9N2	6.0×10^5 TCID ₅₀ /test

Flu A/Yuyao/1/2006	Chicken	H5N1	6.0×10^4 TCID ₅₀ /test
Flu A/Yuyao/2/2006	Chicken	H5N1	6.0×10^4 TCID ₅₀ /test
Flu A/Jiangsu/2/2004	Chicken	H5N1	6.0×10^4 TCID ₅₀ /test
Flu A/Hubei/216/83	Duck	H7N8	3.0×10^5 TCID ₅₀ /test
Flu A/Hubei/118/2003	Duck	H9N2	1.5×10^5 TCID ₅₀ /test
Flu A/Hubei/155/2003	Duck	H9N2	6.0×10^5 TCID ₅₀ /test
Flu A/Hubei/137/1982	Duck	H10N4	3.0×10^4 TCID ₅₀ /test
Flu A/Egypt/3/97	Duck	H5N3	6.0×10^4 TCID ₅₀ /test
Flu A/Henan/1/2004	Tree sparrow	H5N1	6.0×10^5 TCID ₅₀ /test
Flu A/Henan/2/2004	Tree sparrow	H5N1	3.0×10^5 TCID ₅₀ /test
Flu A/Henan/4/2004	Tree sparrow	H5N1	6.0×10^5 TCID ₅₀ /test
Flu A/Wisconsin/66	Turkey	H9N2	6.0×10^4 TCID ₅₀ /test
Flu A/England/1/63	Turkey	H7N3	6.0×10^4 TCID ₅₀ /test
Flu A/Singapore/1/57	Bird	H5N1	6.0×10^4 TCID ₅₀ /test
Flu A/Hunan/71/2004	Bird	H5N1	6.0×10^4 TCID ₅₀ /test
Flu A/S. hanxi/50/2006	Bird	H5N1	6.0×10^4 TCID ₅₀ /test
Flu A/S. hanxi/42/2006	Bird	H5N1	6.0×10^4 TCID ₅₀ /test
Flu A/Fujian/320/2004	Bird	H5N1	3.0×10^5 TCID ₅₀ /test

Influenza A&B Ag Rapid Test can detect all nine influenza B strains.

3. Clinical Study Data Summary

The Influenza A&B Ag Rapid Test performance vs. Cell Culture

Kind of samples	Type	Sensitivity(%)	Specificity(%)	Accuracy(%)
Nasal Swab	A	92.6 (25/27)	96.4 (81/84)	95.5 (106/111)
	B	90.0 (27/30)	95.8 (91/95)	94.4 (118/125)
Throat Swab	A	83.3 (20/24)	95.2 (59/62)	91.9 (79/86)
	B	82.6 (19/23)	91.8 (67/73)	89.6 (86/96)
Nasal Aspirate	A	88.9 (48/54)	93.3 (125/134)	92.0 (173/188)
	B	91.2 (52/57)	95.4 (98/103)	93.8 (150/160)
Nasal discharge/ Nasal mucus	A	80.7 (46/57)	94.9 (93/98)	89.7 (139/155)
	B	89.6 (62/69)	94.6 (87/92)	92.5 (149/161)

4. Analytical Specificity And Cross-reactivity

The Influenza A&B Ag Rapid Test was evaluated with a total of 30 bacterial and viral isolates. Bacterial isolates were evaluated at a concentration between 10^7 and 10^9 org/mL. Viral isolates were evaluated at a concentration of at least 10^4 - 10^8 TCID₅₀/mL. Adenovirus 18 and Parainfluenza virus 3 were tested at 10^2 TCID₅₀/mL. None of the organisms or viruses listed below gave a positive result in the Influenza A&B Ag Rapid Test.

Bacterial Panel:

Acinetobacter calcoaceticus	Bacteroides fragilis
Neisseria gonorrhoeae	Neisseria meningitidis
Pseudomonas aeruginosa	Staphylococcus aureus
Streptococcus pneumoniae	Streptococcus sanguis
Proteus vulgaris	Streptococcus sp. Gp. B
Streptococcus sp. Gp. C	Streptococcus sp. Gp. G
Mycobacterium tuberculosis	Mycoplasma orale

Viral Panel:

Human Adenovirus B	Human Rhinovirus 2
Human Adenovirus C	Human Rhinovirus 14
Adenovirus type 10	Human Rhinovirus 16
Adenovirus type 18	Measles
Human Coronavirus OC43	Mumps
Human Coxsackievirus A9	Sendai virus
Coxsackievirus B5	Parainfluenza virus 2
Human herpesvirus 2	Parainfluenza virus 3

5. Interfering Substances

Whole blood, and several over-the-counter (OTC) products and common chemicals were evaluated and did not interfere with the Influenza A&B Ag Rapid Test at the levels tested: whole blood (2%); three OTC mouthwashes (25%); three OTC throat drops (25%); three OTC nasal sprays (10%); 4-Acetamidophenol (10mg/mL); Acetylsalicylic Acid (20mg/mL); Chlorpheniramine (5mg/mL); Dextromethorphan (10mg/mL); Diphenhydramine (5 mg/mL); Ephedrine (20mg/mL); Guaiacol glyceryl ether (20mg/mL); Oxymetazoline (10mg/mL); Phenylephrine (100 mg/mL); and Phenylpropanolamine (20mg/mL).

For COVID-19 Antigen Rapid Test:

1. Clinical Sensitivity, Specificity and Accuracy

The COVID-19 Antigen Rapid Test Cassette (Swab) has been evaluated with specimens obtained from patients. A commercialized molecular assay was used as the reference method. The results show that the COVID-19 Antigen Rapid Test Cassette (Swab) has a high overall relative accuracy.

Table 1: The COVID-19 Antigen Rapid Test vs PCR

COVID-19 Rapid Test Cassette	Method	PCR		Total Results
	Results	Positive	Negative	
	Positive	39	0	39
	Negative	6	116	122
	Total Results	45	116	161

Relative Sensitivity: 86.7% (95%CI*: 73.2%-95.0%)

*Confidence Intervals

Relative Specificity: 100% (95%CI*: 96.9%-100%)

Accuracy: 96.3% (95%CI*: 92.1%-98.6%)

INDEX OF SYMBOLS

	Consult instructions for use		Tests per kit		Authorized Representative
	For <i>in vitro</i> diagnostic use only		Use by		Do not reuse
	Store between 2-30°C		Lot Number		Catalog#

Göztepe mah. Batışehir cad. No:2/3
İç Kapı No: 1B12 Bağcılar/İSTANBUL
Tel: 00905350499266
www.spesera.com.tr

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
Coronavirus Ag Rapid Test

► COVID-19 Overview

The novel coronaviruses belong to the β genus. COVID-19 is an acute respiratory infectious disease. People are generally susceptible. Currently, patients infected by the novel coronavirus are the main source of infection; asymptomatic infected people can also be an infectious source. The main symptoms include fever, fatigue, loss of smell and dry cough. Nasal congestion, runny nose, sore throat, myalgia and diarrhea are found in a few cases.

SARS viral antigen is generally detectable in upper respiratory specimens during the acute phase of infection. The Coronavirus Ag Rapid Test is a rapid screening tool to detect the presence of SARS viral antigen in the form of a visually interpreted result within minutes.

Product Information

The Coronavirus Ag Rapid Test Cassette (Swab) is an in vitro immunochromatographic assay for the qualitative detection of nucleocapsid protein antigen from SARS-CoV-2 in nasopharyngeal (NP) swab specimens directly or after the swabs has been added to viral transport media from individuals who are suspected of COVID-19 by their healthcare provider. It is intended to aid in the rapid screening of SARS-CoV-2 infections. The Coronavirus Ag Rapid Test Cassette (Swab) does not differentiate between SARS-CoV and SARS-CoV-2.

Specification

Information	InDetail
Time to result	15 minutes
Storage	2-30°C
Shelf life	24 months
Specimen type	Nasopharyngeal swab

Contents

- 25 Test cassettes • 25 Sterile swabs
- 25 Extraction tubes and nozzle
- Workstation • 2 Buffers • Package insert

Test Procedure & Interpretation

Specimen collection


Use the nasopharyngeal swab supplied in the kit.

1. Carefully insert the swab into the nostril of the patient, reaching the surface of posterior nasopharynx. that presents the most secretion under visual inspection.
2. Swab over the surface of the posterior nasopharynx. Rotate the swab several times.
3. Withdraw the swab from the nasal cavity.


Sample preparation

1. Insert the test extraction tube into the workstation in this product. Make sure that the tube is standing firm and reaches the bottom of the workstation.
2. Add 300 µL (about 10 drops) of the extraction buffer into the extraction tube.

3. Insert the swab into the extraction tube which contains 300 µL of the extraction buffer. Roll the swab at least 6 times while pressing the head against the bottom and side of the extraction tube.
4. Leave the swab in the extraction tube for 1 minute.
5. Squeeze the tube several times with fingers from outside of the tube to immerse the swab. Remove the swab. The extracted solution will be used as test sample.
6. Push the nozzle which contains a filter onto the extraction tube. Ensure the nozzle has a tight fit.

Test procedure & interpretation of results

Benefits

- Rapid testing for SARS vial antigen within 15 minutes
- Facilitates patient treatment decisions quickly
- Simple, time-saving procedure
- No instrumentation required
- High sensitivity and specificity

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Sağlık Ürünleri İç ve Dış Tic. Ltd. Şti.