Spesera

Sağlık Ürünleri İç ve Dış Tic. Ltd. Şti.

COVID-19 IgG/IgM Rapid Test

COVID-19 IgG/IgM Rapid Test Cassette (Whole Blood/Serum/Plasma) is a solid phase immunochromatographic assay for the rapid, qualitative and differential detection of IgG and IgM antibodies to 2019 Novel Coronavirus in human whole blood, serum or plasma.



SPESERA SAĞLIK ÜRÜNLERİ İÇ VE DIŞ TİC. LTD. ŞTİ.
Göztepe mah. Batışehir cad. No:2/3
İç Kapı No: 1B12 Bağcılar/İSTANBUL

Tel: 00905350499266
www.spesera.com.tr

Spesera
Sağlık Ürünleri İç ve Dış Tic. Ltd. Şti.

Key Facts of COVID-19

The 2019 new coronavirus pneumonla (NCP), or “COVID-19”, was discoverd for the 2019 Wuhan Viral Pneumonla case In China and was named by world Health Organization on January 12, 2020.

For confirmed coronavirus disease 2019 (COVID-19) cases, reported Illnesses have ranged from mild symptoms to server Illness and death. Symptoms can Include:

- Fever
- Cough
- Shortness of breath


The virus is thought to spread mainly from person-to-person.

- Between peolpe who are in close contact with one another (within about 6 feet).
- Via respiratory droplets produced when an infected person coughs or sneezes.
- These droplets can land in the mouths or noses of people who are nearby or possibly be inhaled into the lungs.

Features & Benefits

- Fast results as soon as 2-10 minutes.
- Facilitates patient treatment decisions quickly.
- Simple, time-saving prodcedure.
- Little specimens, only 5 Ml of serum/plasma or 10 Ml of whole blood specimens.
- All necessary reagents provided & no equipment needed.
- High sensitivity and specificity.

Test Procedure & Interpretation



Ordering Information

Product Description	Specimen	Catalog No.	Format	Kit Size
COVID-19 IgG/IgM Rapid Test	Whole Blood/Serum/Plasma	GCCOV-402a	Cassette	25 Tests/Kit